

INTRODUCTION

The Jes Extender was designed and developed in Denmark. The purpose of the extender is to lengthen the penis by means of traction. The Jes Extender can be used as a gentle mono-therapeutic device. The use of the Jes Extender does not require medical supervision.



The mechanism of the device is based on the principle of traction, a method which simulates the human body's natural reaction to the application of physical force. Tissue cells react by multiplying and the tissue gradually expands, a principle utilized by various indigenous tribes for multiple doubling of the size of earlobes and lips.


The Jes Extender has been tested on patients who have all obtained lengthening of the penis after using the device in accordance with instructions. The lengthening is permanent and does not entail a narrowing of the circumference of the penis.

The test subjects achieved an average lengthening of the erected penis of 2.8 cm (24 per cent) after a treatment period of 3-4 months. The treatment does not affect urination, sexual prowess or the power of reproduction. Lengthening can only be obtained by daily usage of the device and the effect correlates with the number of hours of usage. The above mentioned results were obtained after 12 hours' daily treatment over a period of 3 months.

DESCRIPTION OF THE VARIOUS COMPONENTS AND INSTRUCTION FOR ASSEMBLAGE

The Extender consists of a basal plastic ring (1) with a recess for the urethra (2). The ring is fastened around the root of the penis. Two metal bars, one with a right-handed thread (3), the other with a left-handed thread (4) are attached to the basal ring. Inside each bar, hidden in a tube, is a spring contrivance (5). At the distal end the two bars are attached to a plastic support (6) with two holes, through which the ends of a silicone tube are passed.



The components 1-5 are preassembled on delivery. Depending on the size of the penis either the short or the long extension bars are mounted by means of the round spanners provided (See fig. 3, page 9). The additional extension bars are required for later usage as lengthening of the penis occurs. The plastic support is mounted by introducing the ends of the bars (5) into the two big holes of the support. The support is then pressed towards the basal ring (1). The springs are thus pressed home. Pressure should be exerted until a snap or a click indicates that the support is firmly secured to the bars.



7

INSTRUCTIONS FOR MOUNTING

The device is mounted by fastening the basal ring (1) around the root of the penis, so that the two metal bars (2) run parallel with the penis, one on each side, and the head of the penis rests in the bowl-shaped support (3). By passing the two ends of the silicone tube through the holes in the plastic support (4) a loop is made through which the head of the penis is introduced. Before fastening of the penis the plastic support should be pressed towards the body so that the two springs press home. The silicone tube is fastened around the penis head by pulling the ends of the tube (5) on the reverse side of the support. In men who have not been circumcised the silicone tube is fitted just below the glans with the foreskin in place. The ends are pressed into the two notches (6) to secure the fastening of the penis head. Additional pieces of metal bar (7) may be added as required by means of the two keys (8).



8

PROTECTION PAD ACCESSORY TO THE JES EXTENDER.

In rare cases initial soreness can be experienced from using the Extender. Such soreness is usually caused by initial excess of force, which should not exceed 6-900 gram of traction for 4-8 hours daily during the first 2-3 weeks. The force of the traction should then be increased to 900-1200 grams, and the treatment period to 8-12 hours daily.

To relieve such occasional soreness, the enclosed Protection Pad can be used in the following way:

- Expand the Protection Pad slowly.
- Fit it behind the penis glans, and allow it to contract to a snug fit.
- Mount the extra long silicone tube as usual.

The Protection Pad can be washed in lukewarm water, and new pads can be purchased at your local Extender dealer or from this site.

INSTRUCTIONS FOR ADJUSTMENT

When the device has been fitted the force of traction can be adjusted by turning the two thumb-screws (1) in an inward direction. The actual force exerted can be read on the hexagonal pistons (2) and the force of traction can be varied from 600 to 900, 1200 or 1500 g (see fig. 4). Each marking on the piston equals a traction force of 300 g. Thus, at the exertion of maximum force (1500 g) only one marking is visible on the piston. A force of 600-900 g is recommended for the first 2-3 weeks and during that period the Extender should not be worn for more than 4-5 hours a day. From then on the force may be increased to 900-1200 g and to 12 hours daily usage of the device. The Extender should be worn approximately 12 hours a day for a total period of 3-6 months.

INSTRUCTIONS FOR USE AND DISMOUNTING OF THE EXTENDER

During use of the Extender the penis can be held at any position except sideways. In a pair of loose-fitting trousers the device is inconspicuous. It can be worn at night but may become detached during uneasy sleep or during a nightly erection. Accidental detachment merely requires re-mounting of the device.

The device is dismantled by loosening the silicone tube from the ventral side of the support. The Extender can then be detached. The thumbscrews are screwed home in an outward direction towards the basal ring. When not in use the Extender can be left assembled or it can be disassembled and placed in the box as required.




Fig. 5

1

2

13

MATERIALS USED

1.

Silicone type: WMQ, USP 21-4, PA/HP/exp3/(82)57/1984.

Owing to its chemical composition Silicone rubber is temperature resistant, weather resistant and re-sistant to wear and tear. It is proof against exposure to temperatures ranging from -65°C to $+210^{\circ}\text{C}$.

2 + 4.

Plastic parts: POM (Hostaform), C 2521 G

Proof against temperatures ranging from -20°C to $+80^{\circ}\text{C}$.

3.

Brass: CuZnS9Pb3 - W.no. 2.0401.

Nitin coating: Nitin is an intermetallic composition, combining the properties of nickel and tin, which is produced by electrolytic depositing. Clinical trials have demonstrated that Nitin coatings, despite their contents of 35% nickel, do not seem to precipitate any reaction in nickel allergies.

5.

Stainless steel springs: Din 17224 - W.no. 4310

CLEANING

The JES Extender should be cleaned daily with soap and water.

A chlorhexidine alcohol solution of 0.5% can be used for disinfection.

WORLDWIDE WARRANTY

Your JES Extender is covered by warranty for a period of 12 months under the following conditions: The warranty covers claims which result from defects in material or workmanship. If a defect is discovered the JES Extender must be returned to the dealer together with the completed warranty card.

ENVIROMENTAL ASPECTS

The JES Extender is delivered in packaging indispensable to shipment of the device. Only environment friendly and recyclable materials are used.

Wooden box:

Made of mahogany, harvested from African or North American trees, cut down every ten years as part of routine forestry maintenance.

Glue and varnish:

Water based.

Corrugated paper:

Contains no less than 45% recycled material. Can be recycled.

Plastic bags:

PVC-free. Only carbon dioxide and water are liberated on incineration.

The brochure:

The paper is made of chlorine free bleached cellulose.

Printing ink does not contain heavy metals.

Foam rubber:

Made of PUR foam. Does not emit agents which act on the skin.